

Luciferase Mouse Monoclonal Antibody(7C9)

Catalog No: #44184

Package Size: #44184-1 50ul #44184-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Luciferase Mouse Monoclonal Antibody(7C9)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7C9
Isotype	IgG1
Applications	WB
Species Reactivity	Dm
Specificity	Luciferase Mouse Monoclonal Antibody(7C9) detects endogenous levels of Luciferase
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Luciferase
Other Names	Luciferase
Accession No.	Swiss-Prot:P08659Gene ID:
Uniprot	P08659
Calculated MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:2000

Images



Western blot analysis of Firefly Luciferase Transfected HeLa Cell Lysate using Luciferase Mouse mAb diluted at 1:2000.

Note: This product is for in vitro research use only