

C4b-binding protein

Catalog No: #AP81222



Package Size: #AP81222-1 50ug #AP81222-2 100ug #AP81222-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C4b-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C4bp,C4bp
Accession No.	Uniprot ID:P08607
Uniprot	P08607
GeneID	12269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:C4bpa

Note: This product is for in vitro research use only