

Complement component C1q receptor

Catalog No: #AP81225

Package Size: #AP81225-1 50ug #AP81225-2 100ug #AP81225-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement component C1q receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C1q/MBL/SPA receptor,C1qR(p),Cell surface antigen AA4,Complement component 1 q subcomponent receptor 1,Lymphocyte antigen 68,Ly-68,CD93,Aa4,C1qr1,C1qrp,Ly68
Accession No.	Uniprot ID:O89103
Uniprot	O89103
GeneID	17064;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd93

Note: This product is for in vitro research use only