

Caldesmon

Catalog No: #AP81228



Package Size: #AP81228-1 50ug #AP81228-2 100ug #AP81228-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Caldesmon
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CDM,CAD
Accession No.	Uniprot ID:P12957
Uniprot	P12957
GeneID	373965;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CALD1

Note: This product is for in vitro research use only