

Cadherin-13

Catalog No: #AP81230



Package Size: #AP81230-1 50ug #AP81230-2 100ug #AP81230-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cadherin-13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heart cadherin,H-cadherin,Truncated cadherin,T-cad
Accession No.	Uniprot ID:Q9WTR5
Uniprot	Q9WTR5
GeneID	12554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdh13

Note: This product is for in vitro research use only