

Extracellular calcium-sensing receptor

Catalog No: #AP81236



Package Size: #AP81236-1 50ug #AP81236-2 100ug #AP81236-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Extracellular calcium-sensing receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CaSR, Parathyroid cell calcium-sensing receptor 1, PCaR1, GPRC2A, PCAR1
Accession No.	Uniprot ID:P41180
Uniprot	P41180
GeneID	846;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CASR

Note: This product is for in vitro research use only