

Calcitonin gene-related peptide 1

Catalog No: #AP81237



Package Size: #AP81237-1 50ug #AP81237-2 100ug #AP81237-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Alpha-type CGRP, Calcitonin gene-related peptide I, CGRP-I, CALC, CCALCI
Accession No.	Uniprot ID:Q9MYV1
Uniprot	P41547
GeneID	403946;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CGRP

Note: This product is for in vitro research use only