

Breast cancer type 1 susceptibility protein homolog

Catalog No: #AP81239



Package Size: #AP81239-1 50ug #AP81239-2 100ug #AP81239-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Breast cancer type 1 susceptibility protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RING-type E3 ubiquitin transferase BRCA1
Accession No.	Uniprot ID:O54952
Uniprot	O54952
GeneID	497672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Brca1

Note: This product is for in vitro research use only