

Caspase-4

Catalog No: #AP81245



Package Size: #AP81245-1 50ug #AP81245-2 100ug #AP81245-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CASP-4,ICE and Ced-3 homolog 2,ICH-2,ICE(rel)-II,Mih1,Protease TX,ICH2
Accession No.	Uniprot ID:P49662
Uniprot	P49662
GeneID	837;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CASP4

Note: This product is for in vitro research use only