

Enhancer of filamentation 1

Catalog No: #AP81250

Package Size: #AP81250-1 50ug #AP81250-2 100ug #AP81250-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Enhancer of filamentation 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mEF1,CRK-associated substrate-related protein,CAS-L,Neural precursor cell expressed developmentally down-regulated protein 9,NEDD-9,p105,CasI
Accession No.	Uniprot ID:O35177
Uniprot	O35177
GeneID	18003;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nedd9

Note: This product is for in vitro research use only