

## Cyclin B1 Mouse Monoclonal Antibody(5C1)

Catalog No: #44188

Package Size: #44188-1 50ul #44188-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Cyclin B1 Mouse Monoclonal Antibody(5C1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5C1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Cyclin B1 Mouse Monoclonal Antibody(5C1) detects endogenous levels of CCNB1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Cyclin B1
Other Names	CCNB1
Accession No.	Swiss-Prot:P14635Gene ID:891
Uniprot	P14635
GeneID	891;
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Cyclin B1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Cyclin B1 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only