

Cadherin-1

Catalog No: #AP81260



Package Size: #AP81260-1 50ug #AP81260-2 100ug #AP81260-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cadherin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% by SDS-PAGE
Species Reactivity	Rt
Other Names	Epithelial cadherin,E-cadherin,Uvomorulin,CD324
Accession No.	Uniprot ID:Q9R0T4
Uniprot	Q9R0T4
GeneID	83502;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdh1

Note: This product is for in vitro research use only