

Cadherin-5

Catalog No: #AP81262



Package Size: #AP81262-1 50ug #AP81262-2 100ug #AP81262-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cadherin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Vascular endothelial cadherin,VE-cadherin,CD144
Accession No.	Uniprot ID:Q6URK6
Uniprot	Q6URK6
GeneID	414735;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDH5

Note: This product is for in vitro research use only