

# Cadherin-6

Catalog No: #AP81263



Package Size: #AP81263-1 50ug #AP81263-2 100ug #AP81263-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cadherin-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kidney cadherin,K-cadherin
Accession No.	Uniprot ID:P97326
Uniprot	P97326
GeneID	12563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cdh6

Note: This product is for in vitro research use only