

C4b-binding protein alpha chain

Catalog No: #AP81264

Package Size: #AP81264-1 50ug #AP81264-2 100ug #AP81264-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C4b-binding protein alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	C4bp
Accession No.	Uniprot ID:Q28065
Uniprot	Q28065
GeneID	281651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C4BPA

Note: This product is for in vitro research use only