

E3 ubiquitin-protein ligase CBL-B

Catalog No: #AP81268

Package Size: #AP81268-1 50ug #AP81268-2 100ug #AP81268-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase CBL-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Casitas B-lineage lymphoma proto-oncogene b,RING finger protein 56,RING-type E3 ubiquitin transferase CBL-B,SH3-binding protein CBL-B,Signal transduction protein CBL-B,RNF56,Nbla00127
Accession No.	Uniprot ID:Q13191
Uniprot	Q13191
GeneID	868;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CBLB

Note: This product is for in vitro research use only