

Cyclin-dependent kinase 4

Catalog No: #AP81275

Package Size: #AP81275-1 50ug #AP81275-2 100ug #AP81275-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CRK3, Cell division protein kinase 4, PSK-J3, Crk3
Accession No.	Uniprot ID:P30285
Uniprot	P30285
GeneID	12567;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdk4

Note: This product is for in vitro research use only