

## Cell adhesion molecule 3

Catalog No: #AP81277



Package Size: #AP81277-1 50ug #AP81277-2 100ug #AP81277-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cell adhesion molecule 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Immunoglobulin superfamily member 4B,IgSF4B,Nectin-like protein 1,NECL-1,Igsf4b,Necl1
Accession No.	Uniprot ID:Q1WIM3
Uniprot	Q1WIM3
GeneID	360882;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cadm3

Note: This product is for in vitro research use only