

Alpha-S1-casein

Catalog No: #AP81278



Package Size: #AP81278-1 50ug #AP81278-2 100ug #AP81278-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-S1-casein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-casein,Csn1,Csna
Accession No.	Uniprot ID:P19228
Uniprot	P19228
GeneID	12990;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Csn1s1

Note: This product is for in vitro research use only