

Caspase-12

Catalog No: #AP81283



Package Size: #AP81283-1 50ug #AP81283-2 100ug #AP81283-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CASP-12
Accession No.	Uniprot ID:Q920D5
Uniprot	Q920D5
GeneID	156117;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Casp12

Note: This product is for in vitro research use only