

Corticosteroid-binding globulin

Catalog No: #AP81284



Package Size: #AP81284-1 50ug #AP81284-2 100ug #AP81284-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Corticosteroid-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CBG,Serpin A6,Transcortin,Cbg
Accession No.	Uniprot ID:Q9GK37
Uniprot	Q9GK37
GeneID	396736;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpina6

Note: This product is for in vitro research use only