

## Calpain-2 catalytic subunit

Catalog No: #AP81285



Package Size: #AP81285-1 50ug #AP81285-2 100ug #AP81285-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Calpain-2 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Calcium-activated neutral proteinase 2,CANP 2,Calpain M-type,Calpain large polypeptide L2,Calpain-2 large subunit,Millimolar-calpain,M-calpain,CANPL2
Accession No.	Uniprot ID:P17655
Uniprot	P17655
GeneID	824;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CAPN2

Note: This product is for in vitro research use only