

Transforming growth factor-beta-induced protein ig-h3

Catalog No: #AP81289



Package Size: #AP81289-1 50ug #AP81289-2 100ug #AP81289-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor-beta-induced protein ig-h3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta ig-h3
Accession No.	Uniprot ID:P82198
Uniprot	P82198
GeneID	21810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tgfbi

Note: This product is for in vitro research use only