

Beclin-1 Mouse Monoclonal Antibody(5C5)

Catalog No: #44191

Package Size: #44191-1 50ul #44191-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Beclin-1 Mouse Monoclonal Antibody(5C5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5C5
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Beclin-1 Mouse Monoclonal Antibody(5C5) detects endogenous levels of BECN1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Beclin-1
Other Names	BECN1
Accession No.	Swiss-Prot:Q14457Gene ID:8678
Uniprot	Q14457
GeneID	8678;
Calculated MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Beclin-1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Beclin-1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only