

Caspase-2

Catalog No: #AP81293



Package Size: #AP81293-1 50ug #AP81293-2 100ug #AP81293-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CASP-2, ICH-1L/1S, Protease ICH-1, ICH1
Accession No.	Uniprot ID:Q98943
Uniprot	Q98943
GeneID	395857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CASP2

Note: This product is for in vitro research use only