

Caspase-3

Catalog No: #AP81294



Package Size: #AP81294-1 50ug #AP81294-2 100ug #AP81294-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CASP-3,Apopain,Cysteine protease CPP32, CPP-32,IRP,LICE,Protein Yama,SREBP cleavage activity 1,SCA-1,Cpp32
Accession No.	Uniprot ID:P55213
Uniprot	P55213
GeneID	25402;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Casp3

Note: This product is for in vitro research use only