

## Extracellular calcium-sensing receptor

Catalog No: #AP81296



Package Size: #AP81296-1 50ug #AP81296-2 100ug #AP81296-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Extracellular calcium-sensing receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CaSR, Parathyroid cell calcium-sensing receptor, PCaR1, Gprc2a, Pcar1
Accession No.	Uniprot ID:P48442
Uniprot	P48442
GeneID	24247;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Casr

Note: This product is for in vitro research use only