

Corticosteroid-binding globulin

Catalog No: #AP81297



Package Size: #AP81297-1 50ug #AP81297-2 100ug #AP81297-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Corticosteroid-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CBG,Serpin A6,Transcortin,CBG
Accession No.	Uniprot ID:P08185
Uniprot	P08185
GeneID	866;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPINA6

Note: This product is for in vitro research use only