

## Cyclin-dependent-like kinase 5

Catalog No: #AP81299



Package Size: #AP81299-1 50ug #AP81299-2 100ug #AP81299-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cyclin-dependent-like kinase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cell division protein kinase 5,Serine/threonine-protein kinase PSSALRE,Tau protein kinase II catalytic subunit,TPKII catalytic subunit,Cdkn5
Accession No.	Uniprot ID:Q03114
Uniprot	Q03114
GeneID	140908;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cdk5

Note: This product is for in vitro research use only