

# T-cell surface glycoprotein CD3 delta chain

Catalog No: #AP81305



Package Size: #AP81305-1 50ug #AP81305-2 100ug #AP81305-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	T-cell surface glycoprotein CD3 delta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	T-cell receptor T3 delta chain,CD3d,T3d
Accession No.	Uniprot ID:P04235
Uniprot	P04235
GeneID	12500;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cd3d

Note: This product is for in vitro research use only