

T-cell surface glycoprotein CD3 delta chain

Catalog No: #AP81306



Package Size: #AP81306-1 50ug #AP81306-2 100ug #AP81306-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD3 delta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	T-cell receptor T3 delta chain,CD3d,T3d
Accession No.	Uniprot ID:P19377
Uniprot	P19377
GeneID	25710;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd3d

Note: This product is for in vitro research use only