

C-C chemokine receptor type 7

Catalog No: #AP81310



Package Size: #AP81310-1 50ug #AP81310-2 100ug #AP81310-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	C-C CKR-7,CD197
Accession No.	Uniprot ID:Q5MD62
Uniprot	Q5MD62
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCR7

Note: This product is for in vitro research use only