

C-C motif chemokine 7

Catalog No: #AP81312



Package Size: #AP81312-1 50ug #AP81312-2 100ug #AP81312-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interocrine/chemokine MARC, Monocyte chemoattractant protein 3, Monocyte chemotactic protein 3, MCP-3, Protein FIC, Small-inducible cytokine A7, Fic, Mcp3, Scya7
Accession No.	Uniprot ID:Q03366
Uniprot	Q03366
GeneID	20306;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ccl7

Note: This product is for in vitro research use only