

Cadherin-3

Catalog No: #AP81314



Package Size: #AP81314-1 50ug #AP81314-2 100ug #AP81314-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cadherin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Placental cadherin,P-cadherin,CDHP
Accession No.	Uniprot ID:P22223
Uniprot	P22223
GeneID	1001;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDH3

Note: This product is for in vitro research use only