

Calcitonin gene-related peptide type 1 receptor

Catalog No: #AP81316



Package Size: #AP81316-1 50ug #AP81316-2 100ug #AP81316-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Calcitonin gene-related peptide type 1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CGRP type 1 receptor, Calcitonin receptor-like receptor, Cgrpr
Accession No.	Uniprot ID:Q63118
Uniprot	Q63118
GeneID	25029;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Calcrl

Note: This product is for in vitro research use only