

Cathepsin L2

Catalog No: #AP81318



Package Size: #AP81318-1 50ug #AP81318-2 100ug #AP81318-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin L2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CTSL2
Accession No.	Uniprot ID:Q5E998
Uniprot	Q5E998
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CTSV

Note: This product is for in vitro research use only