

Bone morphogenetic protein receptor type-1A

Catalog No: #AP81320



Package Size: #AP81320-1 50ug #AP81320-2 100ug #AP81320-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein receptor type-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BMP type-1A receptor,Activin receptor-like kinase 3,ALK-3,BMP-2/BMP-4 receptor,Serine/threonine-protein kinase receptor R5,SKR5,CD292,Acvrlk3,Bmpr
Accession No.	Uniprot ID:P36895
Uniprot	P36895
GeneID	12166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Bmpr1a

Note: This product is for in vitro research use only