

Caspase-1

Catalog No: #AP81326



Package Size: #AP81326-1 50ug #AP81326-2 100ug #AP81326-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CASP-1, Interleukin-1 beta convertase, IL-1BC, Interleukin-1 beta-converting enzyme, ICE, IL-1 beta-converting enzyme, p45, IL1BC
Accession No.	Uniprot ID: Q9MZV7
Uniprot	Q9MZV7
GeneID	403724;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: CASP1

Note: This product is for in vitro research use only