

# Cyclin-dependent kinase 12

Catalog No: #AP81330



Package Size: #AP81330-1 50ug #AP81330-2 100ug #AP81330-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cyclin-dependent kinase 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cell division protein kinase 12
Accession No.	Uniprot ID:E1BB50
Uniprot	E1BB50
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CDK12

Note: This product is for in vitro research use only