

Cyclin-dependent kinase 2

Catalog No: #AP81331



Package Size: #AP81331-1 50ug #AP81331-2 100ug #AP81331-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell division protein kinase 2,Cdkn2
Accession No.	Uniprot ID:P97377
Uniprot	P97377
GeneID	12566;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdk2

Note: This product is for in vitro research use only