

Cocaine- and amphetamine-regulated transcript protein

Catalog No: #AP81333

Package Size: #AP81333-1 50ug #AP81333-2 100ug #AP81333-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cocaine- and amphetamine-regulated transcript protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cart
Accession No.	Uniprot ID:P56388
Uniprot	P56388
GeneID	27220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cartpt

Note: This product is for in vitro research use only