

Calretinin

Catalog No: #AP81337



Package Size: #AP81337-1 50ug #AP81337-2 100ug #AP81337-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Calretinin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | CR,29 kDa calbindin,CAB29 |
| Accession No. | Uniprot ID:P22676 |
| Uniprot | P22676 |
| GeneID | 794; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CALB2

Note: This product is for in vitro research use only