

Cadherin-1

Catalog No: #AP81339



Package Size: #AP81339-1 50ug #AP81339-2 100ug #AP81339-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Epithelial cadherin,E-cadherin,Liver cell adhesion molecule,L-CAM
Accession No.	Uniprot ID:P08641
Uniprot	P08641
GeneID	415860;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDH1

Note: This product is for in vitro research use only