

Caprin-1

Catalog No: #AP81346



Package Size: #AP81346-1 50ug #AP81346-2 100ug #AP81346-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Caprin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cytoplasmic activation- and proliferation-associated protein 1, RNA granule protein 105, RNG105
Accession No.	Uniprot ID:Q1LZB6
Uniprot	Q1LZB6
GeneID	535571;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CAPRIN1

Note: This product is for in vitro research use only