

Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein

Catalog No: #AP81351

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP81351-1 50ug #AP81351-2 100ug #AP81351-3 1mg

Description

Product Name	Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-terminal PDZ ligand of neuronal nitric oxide synthase protein, Nitric oxide synthase 1 adaptor protein, Capon, Kiaa0464
Accession No.	Uniprot ID: Q9D3A8
Uniprot	Q9D3A8
GeneID	70729;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nos1ap

Note: This product is for in vitro research use only