

## Dipeptidyl peptidase 1

Catalog No: #AP81355



Package Size: #AP81355-1 50ug #AP81355-2 100ug #AP81355-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dipeptidyl peptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cathepsin C,Cathepsin J,Dipeptidyl peptidase I,DPP-I,DPPI,Dipeptidyl transferase
Accession No.	Uniprot ID:O97578
Uniprot	O97578
GeneID	403458;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CTSC

Note: This product is for in vitro research use only