

# phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4)

Catalog No: #12837

Package Size: #12837-1 50ul #12837-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |  |
|-----------------------|--|
| Product Name          | phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4)                                   |
| Host Species          | Mouse  |
| Clonality             | Monoclonal   |
| Clone No.             | 6B4  |
| Isotype               | IgG1   |
| Applications          | IHC  |
| Species Reactivity    | Hu   |
| Specificity           | phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4) detects endogenous levels of MLKL |
| Immunogen Type        | Peptide  |
| Immunogen Description | Synthetic Peptide  |
| Target Name           | MLKL   |
| Modification          | Phospho  |
| Other Names           | MLKL   |
| Accession No.         | Swiss-Prot:Q8NB16Gene ID:197259  |
| Uniprot               | Q8NB16   |
| GeneID                | 197259;  |
| Calculated MW         | 54kd   |
| Concentration         | 1mg/ml   |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.              |
| Storage               | Store at -20°C/1 year  |


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Phospho-MLKL S358 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using Phospho-MLKL S358 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only