

# phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4)

Catalog No: #12837

Package Size: #12837-1 50ul #12837-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6B4
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	phospho-MLKL (S358) Mouse Monoclonal Antibody(6B4) detects endogenous levels of MLKL
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	MLKL
Modification	Phospho
Other Names	MLKL
Accession No.	Swiss-Prot:Q8NB16Gene ID:197259
Uniprot	Q8NB16
GeneID	197259;
Calculated MW	54kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Phospho-MLKL S358 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using Phospho-MLKL S358 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only