

## E3 ubiquitin-protein ligase CBL-B

Catalog No: #AP81361



Package Size: #AP81361-1 50ug #AP81361-2 100ug #AP81361-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	E3 ubiquitin-protein ligase CBL-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Casitas B-lineage lymphoma proto-oncogene b,RING-type E3 ubiquitin transferase CBL-B,SH3-binding protein CBL-B,Signal transduction protein CBL-B
Accession No.	Uniprot ID:Q8K4S7
Uniprot	Q8K4S7
GeneID	171136;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cblb

Note: This product is for in vitro research use only