

# S-methylmethionine--homocysteine S-methyltransferase BHMT2

Catalog No: #AP81363

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP81363-1 50ug #AP81363-2 100ug #AP81363-3 1mg

## Description

Product Name	S-methylmethionine--homocysteine S-methyltransferase BHMT2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SMM-hcy methyltransferase,Betaine--homocysteine S-methyltransferase 2
Accession No.	Uniprot ID:Q91WS4
Uniprot	Q91WS4
GeneID	64918;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Bhmt2

Note: This product is for in vitro research use only